

Macrophage scavenger receptor types I and II

Catalog No: #AP85444



Package Size: #AP85444-1 50ug #AP85444-2 100ug #AP85444-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage scavenger receptor types I and II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Macrophage acetylated LDL receptor I and II,CD204
Accession No.	Uniprot ID:P21758
Uniprot	P21758
GeneID	281311;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MSR1

Note: This product is for in vitro research use only