

Macrophage scavenger receptor types I and II

Catalog No: #AP85445



Package Size: #AP85445-1 50ug #AP85445-2 100ug #AP85445-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage scavenger receptor types I and II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Macrophage acetylated LDL receptor I and II, Scavenger receptor class A member 1, CD204, SCARA1
Accession No.	Uniprot ID:P21757
Uniprot	P21757
GeneID	4481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MSR1

Note: This product is for in vitro research use only