

Myosin-1

Catalog No: #AP85448



Package Size: #AP85448-1 50ug #AP85448-2 100ug #AP85448-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Myosin heavy chain 1, Myosin heavy chain 2x, MyHC-2x, Myosin heavy chain IIX/d, MyHC-IIX/d, Myosin heavy chain, skeletal muscle, adult 1
Accession No.	Uniprot ID:P12882
Uniprot	P12882
GeneID	4619;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYH1

Note: This product is for in vitro research use only