

NLR family CARD domain-containing protein 4

Catalog No: #AP85451



Package Size: #AP85451-1 50ug #AP85451-2 100ug #AP85451-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NLR family CARD domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Ice protease-activating factor, Ipaf, IPAF
Accession No.	Uniprot ID:F1MHT9
Uniprot	F1MHT9
GeneID	512480;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NLRC4

Note: This product is for in vitro research use only