

Nucleolar GTP-binding protein 1

Catalog No: #AP85453



Package Size: #AP85453-1 50ug #AP85453-2 100ug #AP85453-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleolar GTP-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chronic renal failure gene protein,GTP-binding protein NGB,CRFG,NOG1
Accession No.	Uniprot ID:Q9BZE4
Uniprot	Q9BZE4
GeneID	23560;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GTPBP4

Note: This product is for in vitro research use only