

## Mitogen-activated protein kinase 3

Catalog No: #AP85456



Package Size: #AP85456-1 50ug #AP85456-2 100ug #AP85456-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitogen-activated protein kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase 3,ERT2,Extracellular signal-regulated kinase 1,ERK-1,Insulin-stimulated MAP2 kinase,MAP kinase isoform p44,p44-MAPK,Microtubule-associated protein 2 kinase,p44-ERK1,ERK1,PRKM3
Accession No.	Uniprot ID:P27361
Uniprot	P27361
GeneID	5595;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MAPK3

Note: This product is for in vitro research use only