

Proton myo-inositol cotransporter

Catalog No: #AP85461



Package Size: #AP85461-1 50ug #AP85461-2 100ug #AP85461-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Proton myo-inositol cotransporter |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | H(+)-myo-inositol cotransporter,H(+)-myo-inositol symporter,Solute carrier family 2 member 13 |
| Accession No. | Uniprot ID:Q921A2 |
| Uniprot | Q921A2 |
| GeneID | 171147; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc2a13

Note: This product is for in vitro research use only