

Unconventional myosin-Ib

Catalog No: #AP85463



Package Size: #AP85463-1 50ug #AP85463-2 100ug #AP85463-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Unconventional myosin-Ib
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Myosin I alpha,MMI-alpha,Myosin heavy chain myr 1,Myo1a,Myr1
Accession No.	Uniprot ID:Q05096
Uniprot	Q05096
GeneID	117057;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Myo1b

Note: This product is for in vitro research use only