

# Inositol oxygenase

Catalog No: #AP85464



Package Size: #AP85464-1 50ug #AP85464-2 100ug #AP85464-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inositol oxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Aldehyde reductase-like 6, Kidney-specific protein 32, Myo-inositol oxygenase, MI oxygenase, Renal-specific oxidoreductase, Aldr16, Ksp32, Rsor
Accession No.	Uniprot ID:Q9QXN4
Uniprot	Q9QXN4
GeneID	252899;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene\_full\_name: Miox

Note: This product is for in vitro research use only