

Methylmalonyl-CoA mutase, mitochondrial

Catalog No: #AP85465



Package Size: #AP85465-1 50ug #AP85465-2 100ug #AP85465-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methylmalonyl-CoA mutase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MCM,Methylmalonyl-CoA isomerase
Accession No.	Uniprot ID:Q8MI68
Uniprot	Q8MI68
GeneID	399535;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUT

Note: This product is for in vitro research use only