

Myosin light chain kinase, smooth muscle

Catalog No: #AP85466



Package Size: #AP85466-1 50ug #AP85466-2 100ug #AP85466-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myosin light chain kinase, smooth muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MLCK, Kinase-related protein, KRP, Telokin
Accession No.	Uniprot ID:Q6PDN3
Uniprot	Q6PDN3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mylk

Note: This product is for in vitro research use only