

Kallikrein-10

Catalog No: #AP85469



Package Size: #AP85469-1 50ug #AP85469-2 100ug #AP85469-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kallikrein-10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Normal epithelial cell-specific 1,Protease serine-like 1,NES1,PRSSL1
Accession No.	Uniprot ID:O43240
Uniprot	O43240
GeneID	5655;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KLK10

Note: This product is for in vitro research use only