

# Metallothionein-2A

Catalog No: #AP85475



Package Size: #AP85475-1 50ug #AP85475-2 100ug #AP85475-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metallothionein-2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MT-2A, Metallothionein-IIA, MT-IIA
Accession No.	Uniprot ID:P79379
Uniprot	P79379
GeneID	396827;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MT2A

Note: This product is for in vitro research use only