

Melanoma-derived growth regulatory protein

Catalog No: #AP85479



Package Size: #AP85479-1 50ug #AP85479-2 100ug #AP85479-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanoma-derived growth regulatory protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cartilage-derived retinoic acid-sensitive protein,CD-RAP,Melanoma inhibitory activity protein,Cdrap,Mia1
Accession No.	Uniprot ID:Q62946
Uniprot	Q62946
GeneID	81510;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mia

Note: This product is for in vitro research use only