

# Inositol oxygenase

Catalog No: #AP85482



Package Size: #AP85482-1 50ug #AP85482-2 100ug #AP85482-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inositol oxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Aldehyde reductase-like 6, Myo-inositol oxygenase, MI oxygenase, ALDRL6
Accession No.	Uniprot ID:Q8WN98
Uniprot	Q8WN98
GeneID	397189;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MIOX

Note: This product is for in vitro research use only