

Myelin and lymphocyte protein

Catalog No: #AP85483



Package Size: #AP85483-1 50ug #AP85483-2 100ug #AP85483-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myelin and lymphocyte protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	T-lymphocyte maturation-associated protein,VIP17 proteolipid
Accession No.	Uniprot ID:Q28296
Uniprot	Q28296
GeneID	403933;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAL

Note: This product is for in vitro research use only