

Protein Mdm4

Catalog No: #AP85484



Package Size: #AP85484-1 50ug #AP85484-2 100ug #AP85484-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Protein Mdm4 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Double minute 4 protein,Mdm2-like p53-binding protein,Protein Mdmx,p53-binding protein Mdm4,MDMX |
| Accession No. | Uniprot ID:Q2HJ21 |
| Uniprot | Q2HJ21 |
| GeneID | 514225; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MDM4

Note: This product is for in vitro research use only