

Melanocyte-stimulating hormone receptor

Catalog No: #AP85485



Package Size: #AP85485-1 50ug #AP85485-2 100ug #AP85485-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanocyte-stimulating hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MSH-R,Melanocortin receptor 1,MC1-R,Msh-r
Accession No.	Uniprot ID:Q01727
Uniprot	Q01727
GeneID	17199;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mc1r

Note: This product is for in vitro research use only