

# Mucin-21

Catalog No: #AP85486



Package Size: #AP85486-1 50ug #AP85486-2 100ug #AP85486-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mucin-21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MUC-21,Epiglycanin,C6orf205,UNQ697/PRO1342
Accession No.	Uniprot ID:Q5SSG8
Uniprot	Q5SSG8
GeneID	394263;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MUC21

Note: This product is for in vitro research use only