

Leucyl-cystinyl aminopeptidase

Catalog No: #AP85491



Package Size: #AP85491-1 50ug #AP85491-2 100ug #AP85491-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leucyl-cystinyl aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cystinyl aminopeptidase,Insulin-regulated membrane aminopeptidase,Insulin-responsive aminopeptidase,IRAP,Oxytocinase,OTase,Placental leucine aminopeptidase,P-LAP,OTASE
Accession No.	Uniprot ID:Q9UIQ6
Uniprot	Q9UIQ6
GeneID	4012;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LNPEP

Note: This product is for in vitro research use only