

Stromelysin-2

Catalog No: #AP85492



Package Size: #AP85492-1 50ug #AP85492-2 100ug #AP85492-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stromelysin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SL-2,Matrix metalloproteinase-10,MMP-10,Transin-2
Accession No.	Uniprot ID:O55123
Uniprot	O55123
GeneID	17384;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp10

Note: This product is for in vitro research use only