

Melanin-concentrating hormone receptor 1

Catalog No: #AP85496

Package Size: #AP85496-1 50ug #AP85496-2 100ug #AP85496-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanin-concentrating hormone receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MCH receptor 1,G-protein coupled receptor 24,MCH-1R,MCH1R,SLC-1,Somatostatin receptor-like protein,Gpr24,Slc1
Accession No.	Uniprot ID:P97639
Uniprot	P97639
GeneID	83567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mchr1

Note: This product is for in vitro research use only