

## Mesencephalic astrocyte-derived neurotrophic factor

Catalog No: #AP85499



Package Size: #AP85499-1 50ug #AP85499-2 100ug #AP85499-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mesencephalic astrocyte-derived neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Arginine-rich protein,Protein ARMET,ARMET,ARP
Accession No.	Uniprot ID:P55145
Uniprot	P55145
GeneID	7873;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MANF

Note: This product is for in vitro research use only