

Mitogen-activated protein kinase 14

Catalog No: #AP85501



Package Size: #AP85501-1 50ug #AP85501-2 100ug #AP85501-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	MAP kinase 14,MAPK 14,Mitogen-activated protein kinase p38 alpha,MAP kinase p38 alpha,CSBP1,CSBP2
Accession No.	Uniprot ID:O02812
Uniprot	O02812
GeneID	403856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAPK14

Note: This product is for in vitro research use only