

Induced myeloid leukemia cell differentiation protein Mcl-1

Catalog No: #AP85504

Package Size: #AP85504-1 50ug #AP85504-2 100ug #AP85504-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Induced myeloid leukemia cell differentiation protein Mcl-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Bcl-2-like protein 3,Bcl2-L-3,Bcl-2-related protein EAT/mcl1,mcl1/EAT,BCL2L3
Accession No.	Uniprot ID:Q07820
Uniprot	Q07820
GeneID	4170;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MCL1

Note: This product is for in vitro research use only