

Lactoylglutathione lyase

Catalog No: #AP85507



Package Size: #AP85507-1 50ug #AP85507-2 100ug #AP85507-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lactoylglutathione lyase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Aldoketomutase, Glyoxalase I, Glx I, Ketone-aldehyde mutase, Methylglyoxalase, S-D-lactoylglutathione methylglyoxal lyase
Accession No.	Uniprot ID:Q6P7Q4
Uniprot	Q6P7Q4
GeneID	294320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Glo1

Note: This product is for in vitro research use only