

Microphthalmia-associated transcription factor

Catalog No: #AP85510



Package Size: #AP85510-1 50ug #AP85510-2 100ug #AP85510-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Microphthalmia-associated transcription factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class E basic helix-loop-helix protein 32,bHLHe32,BHLHE32
Accession No.	Uniprot ID:O75030
Uniprot	O75030
GeneID	4286;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MITF

Note: This product is for in vitro research use only