

Myotubularin

Catalog No: #AP85511



Package Size: #AP85511-1 50ug #AP85511-2 100ug #AP85511-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myotubularin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Phosphatidylinositol-3,5-bisphosphate 3-phosphatase,Phosphatidylinositol-3-phosphate phosphatase
Accession No.	Uniprot ID:Q6AXQ4
Uniprot	Q6AXQ4
GeneID	288762;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mtm1

Note: This product is for in vitro research use only