

# Muscle, skeletal receptor tyrosine-protein kinase

Catalog No: #AP85512



Package Size: #AP85512-1 50ug #AP85512-2 100ug #AP85512-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Muscle, skeletal receptor tyrosine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Muscle-specific tyrosine-protein kinase receptor,MuSK
Accession No.	Uniprot ID:O15146
Uniprot	O15146
GeneID	4593;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MUSK

Note: This product is for in vitro research use only