

Protein NDNF

Catalog No: #AP85520



Package Size: #AP85520-1 50ug #AP85520-2 100ug #AP85520-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein NDNF
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neuron-derived neurotrophic factor,C4orf31,UNQ2748/PRO6487
Accession No.	Uniprot ID:Q8TB73
Uniprot	Q8TB73
GeneID	79625;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NDNF

Note: This product is for in vitro research use only