

# Lysosome-associated membrane glycoprotein 1

Catalog No: #AP85522

Package Size: #AP85522-1 50ug #AP85522-2 100ug #AP85522-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysosome-associated membrane glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LAMP-1,120 kDa lysosomal membrane glycoprotein,LGP-120,CD107 antigen-like family member A,CD107a,Lamp-1
Accession No.	Uniprot ID:P14562
Uniprot	P14562
GeneID	25328;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lamp1

Note: This product is for in vitro research use only