

L-lactate dehydrogenase A chain

Catalog No: #AP85523



Package Size: #AP85523-1 50ug #AP85523-2 100ug #AP85523-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase A chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	LDH-A,LDH muscle subunit,LDH-M
Accession No.	Uniprot ID:P00339
Uniprot	P00339
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LDHA

Note: This product is for in vitro research use only