

# NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial

Catalog No: #AP85524

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP85524-1 50ug #AP85524-2 100ug #AP85524-3 1mg

## Description

Product Name	NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Complex I-20kD, CI-20kD, NADH-ubiquinone oxidoreductase 20 kDa subunit, PSST subunit
Accession No.	Uniprot ID:P42026
Uniprot	P42026
GeneID	338079;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: NDUFS7

Note: This product is for in vitro research use only