

Neurabin-2

Catalog No: #AP85525



Package Size: #AP85525-1 50ug #AP85525-2 100ug #AP85525-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurabin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neurabin-II, Neural tissue-specific F-actin-binding protein II, PP1bp134, Protein phosphatase 1 regulatory subunit 9B, Spinophilin, p130
Accession No.	Uniprot ID:O35274
Uniprot	O35274
GeneID	84686;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ppp1r9b

Note: This product is for in vitro research use only