

Neuromodulin

Catalog No: #AP85528



Package Size: #AP85528-1 50ug #AP85528-2 100ug #AP85528-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuromodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Axonal membrane protein GAP-43, Growth-associated protein 43, Protein F1
Accession No.	Uniprot ID:P07936
Uniprot	P07936
GeneID	29423;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gap43

Note: This product is for in vitro research use only