

Nuclear factor erythroid 2-related factor 1

Catalog No: #AP85530



Package Size: #AP85530-1 50ug #AP85530-2 100ug #AP85530-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor erythroid 2-related factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NF-E2-related factor 1,Locus control region-factor 1,Nuclear factor, erythroid derived 2, like 1,Transcription factor 11,TCF-11,Transcription factor HBZ17,Transcription factor LCR-F1,HBZ17,NRF1,TCF11
Accession No.	Uniprot ID:Q14494
Uniprot	Q14494
GeneID	4779;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFE2L1

Note: This product is for in vitro research use only