

Membrane cofactor protein

Catalog No: #AP85531



Package Size: #AP85531-1 50ug #AP85531-2 100ug #AP85531-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Membrane cofactor protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD46,MCP
Accession No.	Uniprot ID:Q6VE48
Uniprot	Q6VE48
GeneID	280851;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD46

Note: This product is for in vitro research use only