

Mitogen-activated protein kinase 14

Catalog No: #AP85533



Package Size: #AP85533-1 50ug #AP85533-2 100ug #AP85533-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase 14,CRK1,Mitogen-activated protein kinase p38 alpha,MAP kinase p38 alpha,Crk1,Csbp1,Csbp2
Accession No.	Uniprot ID:P47811
Uniprot	P47811
GeneID	26416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mapk14

Note: This product is for in vitro research use only