

Membrane metallo-endopeptidase-like 1

Catalog No: #AP85534



Package Size: #AP85534-1 50ug #AP85534-2 100ug #AP85534-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Membrane metallo-endopeptidase-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NEP2(m),Neprilysin II,NEPII,Neprilysin-2,NEP2,Mell1,Nep2
Accession No.	Uniprot ID:P0C1T0
Uniprot	P0C1T0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmel1

Note: This product is for in vitro research use only