

Myosin-11

Catalog No: #AP85537



Package Size: #AP85537-1 50ug #AP85537-2 100ug #AP85537-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myosin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Myosin heavy chain 11, Myosin heavy chain, smooth muscle isoform, SMMHC
Accession No.	Uniprot ID:Q63862
Uniprot	Q63862
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Myh11

Note: This product is for in vitro research use only