

Myelin proteolipid protein

Catalog No: #AP85538



Package Size: #AP85538-1 50ug #AP85538-2 100ug #AP85538-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myelin proteolipid protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PLP,Lipophilin,PLP
Accession No.	Uniprot ID:Q712P7
Uniprot	Q712P7
GeneID	397029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLP1

Note: This product is for in vitro research use only