

# Neurogranin

Catalog No: #AP85539



Package Size: #AP85539-1 50ug #AP85539-2 100ug #AP85539-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurogranin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ng,Protein kinase C substrate 7.5 kDa protein,RC3
Accession No.	Uniprot ID:Q04940
Uniprot	Q04940
GeneID	64356;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nrgn

Note: This product is for in vitro research use only