

## Nuclear factor interleukin-3-regulated protein

Catalog No: #AP85540



Package Size: #AP85540-1 50ug #AP85540-2 100ug #AP85540-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear factor interleukin-3-regulated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	E4 promoter-binding protein 4, Embryo implantation-related NFIL3/E4BP4-like transcription factor, E4bp4
Accession No.	Uniprot ID:O08750
Uniprot	O08750
GeneID	18030;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Nfil3

Note: This product is for in vitro research use only