

Protein NKG7

Catalog No: #AP85541



Package Size: #AP85541-1 50ug #AP85541-2 100ug #AP85541-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein NKG7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Natural killer cell protein 7
Accession No.	Uniprot ID:Q2KJ11
Uniprot	Q2KJ11
GeneID	513006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NKG7

Note: This product is for in vitro research use only