

Mitogen-activated protein kinase kinase kinase kinase 1

Catalog No: #AP85544

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85544-1 50ug #AP85544-2 100ug #AP85544-3 1mg

Description

Product Name	Mitogen-activated protein kinase kinase kinase kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hematopoietic progenitor kinase,MAPK/ERK kinase kinase kinase 1,MEK kinase kinase 1,HPK1
Accession No.	Uniprot ID:Q92918
Uniprot	Q92918
GeneID	11184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAP4K1

Note: This product is for in vitro research use only