

Cell surface glycoprotein MUC18

Catalog No: #AP85545



Package Size: #AP85545-1 50ug #AP85545-2 100ug #AP85545-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cell surface glycoprotein MUC18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gicerin,Melanoma cell adhesion molecule,Melanoma-associated antigen MUC18,CD146,Muc18
Accession No.	Uniprot ID:Q8R2Y2
Uniprot	Q8R2Y2
GeneID	84004;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mcam

Note: This product is for in vitro research use only