

## NKG2D ligand 2

Catalog No: #AP85547

Package Size: #AP85547-1 50ug #AP85547-2 100ug #AP85547-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NKG2D ligand 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N2DL-2,ALCAN-alpha,Retinoic acid early transcript 1H,UL16-binding protein 2,N2DL2,RAET1H,UNQ463/PRO791
Accession No.	Uniprot ID:Q9BZM5
Uniprot	Q9BZM5
GeneID	80328;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ULBP2

Note: This product is for in vitro research use only