

Pyruvate kinase PKLR

Catalog No: #AP85549



Package Size: #AP85549-1 50ug #AP85549-2 100ug #AP85549-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyruvate kinase PKLR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	L-PK,Pyruvate kinase isozymes L/R
Accession No.	Uniprot ID:P12928
Uniprot	P12928
GeneID	24651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PkI

Note: This product is for in vitro research use only