

NEDD8 ultimate buster 1

Catalog No: #AP85555



Package Size: #AP85555-1 50ug #AP85555-2 100ug #AP85555-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NEDD8 ultimate buster 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Negative regulator of ubiquitin-like proteins 1, Renal carcinoma antigen NY-REN-18, NYREN18
Accession No.	Uniprot ID:Q9Y5A7
Uniprot	Q9Y5A7
GeneID	51667;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NUB1

Note: This product is for in vitro research use only