

Collagenase 3

Catalog No: #AP85564



Package Size: #AP85564-1 50ug #AP85564-2 100ug #AP85564-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collagenase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Matrix metalloproteinase-13,MMP-13
Accession No.	Uniprot ID:P33435
Uniprot	P33435
GeneID	17386;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp13

Note: This product is for in vitro research use only